



**PT-134/55**

**VACUUM FILTER HOLDERS MANIFOLD**



**PT-134-05**

**VACUUM FILTER HOLDERS**



Cat. No.	Item's Description	Rate
PT-134	These mechanism are used in microbiological quality control e.g. of water, beer, wine , soft drinks & other laboratory reagents. Analytical determination e.g. of metal particles in lubricating oil & iron oxides in boiler water. Stainless steel and borosilicate glass frit filter support gives optionally uniform distribution of retained bacteria or particles on the memberane filter surface. It can be flamed to sanitize between samples without membranes.	
PT-134/05	<b>STAINLESS STEEL VACUUM FILTER HOLDER</b> consists of lid with teflon sealing ring, 250ml funnel base clamp to connect filter section, detachable 47mm filter disk, silicone gasket ring with rubber cork to connect with 1000 ml. borosilicate filter flask.	<b>5950-00 each</b>
PT-134/10	<b>VACUUM FILTER HOLDER, SYRINGE PATTERN</b> for filtering & cleaning small volumes of liquids, 25mm disk without membrane.	<b>975-00 each</b>
PT- 134/15	<b>BOROSILICATE GLASS VACUUM FILTER HOLDER</b> consisting of funnel reservoir graduated in increments fixed sintered disk in filtration section with special grinding to support membrane, anodized aluminum clamp & heavy wall filter flask.	

	<b>DISK DIA MM</b>	<b>FUNNEL CAP ML</b>	<b>FILTER FALSK CAP ML</b>	<b>STANDARD JOINT/ RUBBER FITTING</b>	<b>RATE</b>
PT-134/20	25	25	250	Rubber	<b>1800-00 each</b>
PT-134/25	25	50	500	Rubber	<b>1900-00 each</b>
PT-134/30	47	300	1000	Rubber	<b>3375-00 each</b>
PT-134/35	25	25	250	Standard	<b>1875-00 each</b>
PT-134/40	25	50	500	Standard	<b>1975-00 each</b>
PT-134/45	47	300	1000	Standard	<b>3525-00 each</b>
PT-134/50	47	300	2000	Standard	<b>4075-00 each</b>