

PT-1050	DROP POINT OF GREASE APPARATUS (As per IP 31, IP 132 and ASTM D-566 specifications) As per IP 31 the apparatus consists of brass sleeve and case, metal cup, a glass boiling tube, cork fitted to a bath and a 600 ml. beaker. The assembly is kept on a heater with energy regulator control, stand and a F.H.P. electric stirrer. As per IP 132 the unit comprises of a glass beaker with fibre cover, glass test tube with three indentations, chromium plated brass cup, cork for tube and a 1.5 mm. diameter brass rod. The assembly is kept on a heater with energy regulator, stand and a F.H.P. electric stirrer.	13500-00 each
PT-1052	Same as PT-1050 but As per IP 132 & ASTM D-566	13500-00 each
PT-1052/01	Same as above with electronic digital temp. Indicator & Controller	19500-00 each

OPTIONAL ACCESSORIES

PT-1052/05	IP 40C low range thermometer.	475-00 each
PT-1052/10	IP 41C high range thermometer.	475-00 each
PT-1052/15	IP 62C drop point special thermometer.	475-00 each