



PT - 11 THIN LAYER CHROMATOGRAPHY APPARATUS

Cat. No.	Item's Description	Rate
PT-11	Thin Layer Chromatography Apparatus (Applicator) made of brass for 0-2mm thickness, Perspex base, Plate rack aluminium, developing tank, Spotting template, T.L.C Plates set, Glass sprayer set, Micro-pipette, Scribe.	9975-00- each
PT-11A	Thin Layer Chromatography App. aluminum applicator. ACCESSORIES FOR T.L.C. APPARATUS	8950-00 each
PT-11/01	BRASS APPLICATOR Spare	4500-00 each
PT-11/02	ALUMINUM APPLICATOR Spare	3200-00 each
PT-11/05	SPOTTING TEMPLATE Made of perspex.	525-00 each
PT-11/10	PLATE RACK Made of anodised aluminum to take ten plates.	980-00 each
PT-11/25	GLASS PLATES All sides ground, for T.L.C 20 x 20cm	55-00 each
PT-11/26	----- do ----- 20 x 10 cm	38-00 each
PT-11/27	----- do ----- 20 x 5 cm	30-00 each
PT-11/30	CHROMATOGRAPHY TANK With cover for plates of different sizes. 20 x 10cm	2250-00 each
PT-11/31	----- do ----- 20 x 20 cm	2500-00 each