



## PT- 1155 E DISINTEGRATION TEST APPARATUS

Cat. No.	Item's Description	Rate
PT-1155	<p><b>DISINTEGRATION TEST APPARATUS</b> I.P. Standard, USP and B.P. Standard.</p> <p>The instrument is designed for accurate estimation of disintegration time of tablets/capules. It is assembled in painted sheet metal cabinet with chemical resistant powder coating, fitted with superior quality constantspeed sturdy motor. The guide discs &amp; strokes are confirming to pharmacopoeial standards.</p> <p><b>Specification :</b></p> <p>Displacement : 50 to 60 mm</p> <p>RPM : 28 to 32 RPM</p> <p>Motor Single : 1/8 HP</p> <p>Heater Single : 125 Watt</p>	<b>14950-00 each</b>
PT- 1155 A	<b>DISINTEGRATION TEST APPARATUS</b> Provided with double basket system : 1/6 H. P. Motor and 250 watts Heater.	<b>18950-00 each</b>
PT- 1155 B	<b>DISINTEGRATION TEST APPARATUS</b> Provided with digital timer and time hold facility in single unit	<b>20950-00 each</b>
PT- 1155 C	<b>DISINTEGRATION TEST APPATATUS</b> Provided with digital timer and time hold facility in double unit	<b>26950-00 each</b>
PT- 1155 D	<b>DISINTEGRATION TEST APPATATUS</b> Provided With Lcd Display Complet with double display of temperature, Auto Parking of Basket, Time Totolicer and lower water level cut of heater in double unit	<b>49750-00 each</b>
PT- 1155 E	<b>DISINTEGRATION TEST APPARATUS</b> same as PT-1155 D but provided with 304 grade stainless steel cabinet (GMP Model)	<b>59750-00 each</b>
PT- 1155 F	<b>SPARE BASKET WITH DISK &amp; GLASS TUBES</b>	<b>2150-00 each</b>
PT- 1155 G	<b>SPARE GLASS TUBE</b>	<b>90-00 each</b>
PT- 1155 H	<b>SPARE DISK</b>	<b>60-00 each</b>