



PT-275 SEROLOGICAL WATER BATH

Cat. No.	Item's Description	Rate
PT-275	<p>SEROLOGICAL WATER BATH Serological water bath is suitable for clinical procedure, incubation, inactivation & pharmaceutical biomedical procedure.</p> <p>CONSTRUCTIONS : Double walled chambers made of stainless steel & outer wall made of thick steel finished with epoxy based powder coated paint. The double walled chambers is provided with a drain plug to facilitate easy emptying & cleaning of the inner chamber when required S.S. Lid & perforated removable diffuser are also provided with double walled chamber. This chamber house an ISI mark 1500 watts heating device, which is controlled by capillary type thermostat from ambient to 100°C within $\pm 1^\circ\text{C}$, two indicator lights & on/off Rotary Switch.</p>	
SIZES OF INNER CHAMBERS		
PT-275/05	300x250x175mm to accommodate	2 racks. 7500-00 each
PT-275/10	350x300x175mm to accommodate	4 racks. 9950-00 each
PT-275/15	450x300x175mm to accommodate	6 racks. 11500-00 each
PT-275-A	SEROLOGICAL WATER BATH provided with electronic digital temperature indicator cum controller.	
PT-275A/05	300x250x175mm to accommodate	2 racks. 15950-00 each
PT-275A/10	350x300x175mm to accommodate	4 racks. 18450-00 each
PT-275A/15	450x300x175mm to accommodate	6 racks. 19950-00 each
PT-275-B	SEROLOGICAL WATER BATH provided with 304 grade Stainless Steel outer body with digital temperature indicator cum controller.	
PT-275B/05	300x250x175mm to accommodate	2 racks. 19950-00 each
PT-275B/10	350x300x175mm to accommodate	4 racks. 22450-00 each
PT-275B/15	450x300x175mm to accommodate	6 racks. 24950-00 each
PT-275-C	SEROLOGICAL WATER BATH fitted with 1/20 H. P. motor with Stainless Steel rod & blades at an extra cost.	4950-00 each