

PT-457

**AUTO CLAVE VERTICAL ELECTRICALLY** heated vertical autoclave is used for sterilization of surgical instruments, dressing material, linen, rubber plastic material by means of saturated steam under pressure of 15 to 20 PSI. It is also used in scientific research institutions & universities.

**CONSTRUCTION:** Sturdy double walled construction with inner chamber is made of 304 grade Stainless Steel 16 SWG sheet & outer shell is also made of Stainless Steel highly polished. The boiler & outer shell is provided with air insulation. Lid is made of thick gauge Stainless Steel sheet & is tightened all around by wing nuts & jointless Silicon gasket. Autoclave is fitted with water level arrangement, dial pressure gauge, steam release valve, spring loaded safety valve; which can be set at any selected point from 10 psi to 20 psi  $\pm$  3 psi, dead weight type safety valve, ISI marked immersion type heating element & Stainless Steel basket with stand.

	<i>Size mm.</i>	<i>Capacity Ltr.</i>	<i>Rate</i>
PT-457/05	300x550	40	<b>27500-00 each</b>
PT-457/10	350x550	50	<b>34500-00 each</b>
PT-457/15	400x600	100	<b>64500-00 each</b>