

PT-480	RESEARCH POLARIMETER Used in medical, pharmaceutical and chemical research laboratories. Consists of a imported half shade nicol polariser, analyser, 100 mm. diameter circular scale divided in degrees 0-360° x 1° and a sugar scale divided –30° to + 130° x 1° ISS, a divided micrometer drum is given to read 0.05° (3 minutes) and 0.1° ISS, the scale is properly illuminated for taking readings, the body is very sturdy and separate sodium lamp stand is provided. Complete with polarimeter tube, bubble traps at the centre, 200 mm. sodium lamp with stand and transformer. 29500-00 each
PT-485	Spare sodium lamp, and transformer for above. 5950-00 each
PT-490	RESEARCH POLARIMETER Same as PT-480 but with 400 mm. tube with sodium lamp, stand and transformer. 39950-00 each

POLARIMETER TUBES

PT-500	Screw caps and cover glasses. Packed in velvet case.	
PT-500/05	Length 100 mm. with centre bulb	925-00 each
PT-500/10	Length 200 mm. with centre bulb	950-00 each
PT-500/20	Length 400 mm. with centre bulb	1050-00 each
PT-505	POLARIMETER TUBES With centre opening, with metallic screw caps and cover glasses. Packed in velvet case.	
PT-505/05	Length 100 mm. with centre cup	950-00 each
PT-505/10	Length 200 mm. with centre cup	975-00 each
PT-505/20	Length 400 mm. with centre cup	1075-00 each