

# Rotary Vacuum Evaporator



## RAPIDVAP SERIES

VERTICAL CONDENSER

**MODEL: PT-VCT**

**PT-70**

Cat. No.	Item's Description	Rate
<b>PT-70</b>	This model of evaporator is supplied complete with geared drive non sparking noiseless motor offering smooth speed control by analog knob upto 200RPM. Jack arrangement is fitted on tapered front base for positioning glass assembly. Double walled <b>water bath thermostatically controlled temp</b> , <b>Vertical</b> glass condenser double coiled with three opening for vacuum and cooling. pear (kjeldhal) shaped evaporating flask capacity 800 ml, receiver flask capacity one litre fitted with a ball joint, feeding tube, vapour path tube provided with a seal, clips, & spring a <b>complete working unit</b> .	<b>Rs. 51950-00</b>

### Replacement parts for PT-70

<b>PT-70/05</b>	Vertical condenser double coiled with two outlets for water circulation. With S-35 Ball at bottom for connecting to receiving flask. quickfit type joint at top for vacuum. cock for feeding tube on one side while the other side will be connected to motor with ring.	<b>11000-00</b>
PT-70/10	Cone Adapter for condenser to Vacuum pump	<b>695-00</b>
PT-70/15	Vapour path tube B-24. B-29	<b>895-00</b>
PT-70/20	Gasket seal for vapour tube (VITON/TEFLON)	<b>695-00</b>
PT-70/25	Gasket Seal for specific Chemicals	<b>300-00</b>
PT-70/30	Feed tube for condenser without Teflon pipe	<b>295-00</b>
PT-70/35	Receiving flask with S-35 cup (socket) joint capacity 1000ml.	<b>850-00</b>
PT-70/40	Evaporating flask Kjeldhal Cap. 800 ml. B-24	<b>850-00</b>
PT-70/45	Evaporating flask Kjeldhal Cap. 500 ml. B-24	<b>695-00</b>
PT-70/50	Evaporating flask Kjeldhal Cap. 300 ml. B-24	<b>595-00</b>
Pt-70/55	Evaporating flask Kjeldhal Cap. 100 ml. B-24	<b>495-00</b>
PT-70/60	Tension spring to hold condensor	<b>90-00</b>
PT-70/65	Wire Clip to hold vapour tube and boiling flask	<b>55-00</b>
PT-70/70	Ball Joint Clip S-35 with screw lock	<b>325-00</b>
PT-70/75	Waterbath double walled inner SS With thermostatic temperature control	<b>7500-00</b>
PT-70/80	Spare heater for water bath (ISI Mark)	<b>675-00</b>
PT-70/85	Oil free vacuum pump with vacuum gauge & regulator	<b>19500-00</b>
PT-70/90	Oil free / chemical resistant vacuum pump with vacuum gauge & Regulator	<b>24950-00</b>
PT-70/95	Glass safety trap	<b>1395-00</b>

Rotary Evaporator - Scientist's first choice since 1960

