

# Rotary Vacuum Evaporator



## RAPIDVAP SERIES

COLD TRAP CONDENSER

**MODEL: PT-CTCT**

**PT-72**

Cat. No.	Item's Description	Rate
<b>PT-72</b>	This model of evaporator is supplied complete with geared drive non sparking noiseless motor offering smooth speed control by analog knob up to 200RPM . Jack arrangement is fitted on tapered front base for positioning glass assembly. Double walled <b>water bath with Thermostatic temperature controller</b> ,glass condenser double walled with one opening for vacuum & Top reservoir for cooling medium, pear (kjeldhal) shaped evaporating flask capacity 800 ml, receiver flask capacity one litre fitted with a ball joint, feeding tube, vapour path tube provided with a seal, spring & clips <b>a complete working unit.</b>	<b>Rs. 52950-00</b>

### Replacement parts for PT-72

PT-72/05	Cold Trap condenser double walled with S-35 ball at bottom for connecting to receiving flask, with one opening for vacuum & top reservoir for cooling medium cock for feeding on one side while the other side will be connected to motor with ring.	<b>12000-00</b>
PT-72/10	Vapour path tube B-24 B-29	<b>895-00</b>
PT-72/15	Gasket seal for vapour tube (Viton/Teflon)	<b>695-00</b>
PT-72/20	Gasket seal for Specific chemicals	<b>300-00</b>
PT-72/25	Feed tube for condenser without Teflon pipe	<b>295-00</b>
PT-72/30	Receiving flask with S-35 cup socket Joint Capacity 1000ml	<b>850-00</b>
PT-72/35	Evaporating Flask Kjeldhal Cap.800ml-B-24	<b>850-00</b>
PT-72/40	Evaporating Flask Kjeldhal Cap.500ml-B-24	<b>695-00</b>
PT-72/45	Evaporating Flask Kjeldhal Cap.300ml-B-24	<b>595-00</b>
PT-72/50	Evaporating Flask Kjeldhal Cap.100ml-B-24	<b>495-00</b>
PT-72/55	Tension spring to hold condenser	<b>90-00</b>
PT-72/60	Wire Clip	<b>55-00</b>
PT-72/65	Ball Joints Clip S-35 with Screw lock	<b>325-00</b>
PT-72/70	Water bath double walled inner SS with thermostatic temperature controller	<b>7500-00</b>
PT-72/75	Spare Heater for water Bath (ISI Mark)	<b>675-00</b>
PT-72/80	Oil free vacuum pump with vacuum gauge & Regulator	<b>19500-00</b>
<b>PT-72/85</b>	<b>Oil free / chemical resistant vacuum pump with vacuum gauge &amp; regulator</b>	<b>24950-00</b>
PT-72/90	Glass Safety Trap	<b>1395-00</b>

Rotary Evaporator - Scientist's first choice since 1960

