



PT-80
CHILLER
REFRIGERATED CIRCULATOR

Model : PT-6.5L

PT-80 CHILLER Refrigerated Circulator has built in circulating pump and water tank of rust proof stainless steel. The built in pump has a flow rate of 2.5 meters and internal tank capacity 6.5 liters. The temperature is controlled by a digital temp. indicator cum controller.

Construction The body of the chiller is mounted on a sturdy angle frame with inner chamber consisting of thick polished stainless steel reservoir with cooling coil to ensure perfect cooling from -10°C to below ambient : The hollow space between the inner chamber and outer wall is filled with high grade puff insulation to prevent loss of temperature. The refrigeration system consisting of hermetically sealed ISI marked reputed make compressor for minimum noise level and smooth operation.

Rs. 69950-00