

PT-800	<b>pH METER</b> With built-in electronic buffer.	
	Readout	: 3½ digit LED display
	Range	: pH: 0 to 14.00                      mV: 0 to ± 1.999
	Resolution	: pH: 0.01pH                              mV: 1
	Accuracy	: pH: ± 0.02 pH ± 2 least count mV: ± 0.2% + 1 least count
	Temp. comp.	: 0 to 100° C
	Power required	: 230V AC ± 10% 50 Hz, 4 VA
	Warming time	: 5 minutes
	Weight	: 1Kg (approx.)                              Body: ABS
	Dimension	: 80 mm (H) x 235 mm (W) x 150 mm (L)
	Electrode	: Combined PVC sleeved
	Electronic buffer	: Equivalent signal output of ideal electrode in solution of 9.2, 7.00, 4.00 pH. Useful for self diagnostic.
	Accessories	: Stand set, Patch cord for Electronic buffer, Buffer tab lets for 9.20, 7.00, 4.00 pH. Screwdriver and Dust proof cover.
PT-800/05	Digital pH Meter (0-14pH) with UNBREAKABLE COMBINED ELECTRODE & Built in Electronic Buffer.	<b>14950-00 each</b>
PT-800/10	Digital pH Meter (0-14pH) with UNBREAKABLE COMBINED ELECTRODE & Thermometer, Auto Temperature Compensation; Electronic Buffer with Temperature Probe.	<b>17950-00 each</b>
PT-810	Spare P.H. Combined Electrode 0-14 PH	<b>975-00 each</b>
PT-820	Stand for electrode	<b>375-00 each</b>