



PT- 859

DIGITAL COLONY COUNTER

Cat. No.	Item's Description	Rate
PT-859	<p>DIGITAL COLONY COUNTER : Colony counter is ideal for rapid , accurate counting and recording of bacterial & mould colonies in biology, zoology bacteriology, pathology laboratories of any place where colony counting is a regular routine.</p> <p>Colony countergrid is illuminated with fluorescent lamp for uniform cool light illumination and is provided with a magnifying glass of ease of viewing special pentype probe allows marking on petridish with an audio beep and counting simultaneously. Counting is registered on 4 digit resettable digital counter.</p> <p>Counter : 4 digit (9999 counts) resettable digital counter with audio beep at every count.</p> <p>Probe : Auto Marker Pen probe for marking and counting simultaneously.</p> <p>Grid : Wolffuegel Grid glass plate.</p> <p>Illumination : Uniform glare free illumination by cool fluorescent tube light.</p> <p>Magnification : 100mm. diameter Magnifying Lens.</p> <p>Dimension : 250 W x 310 D x 13 mm H.</p> <p>Weight : 4.5 kg Approx</p>	9950-00 each
PT-860	PETRI PLATTER For above (Manual).	675-00 each