



**PT-250**  
**WATER BATH THERMOSTATIC**

Cat. No.	Item's Description	Rate
<b>WATER BATHS</b>		
PT-250	<b>WATER BATH THERMOSTATIC</b> Rectangular water bath is suitable for evaporating liquids like alcohol or for melting solids such as agar. Concentric rings can be removed to accommodate containers of various sizes.	
	<b>CONSTRUCTIONS</b> Double walled chambers made of stainless steel and outer wall made of thick mild steel finished with epoxy based powder coated paint. The double walled chamber is filled with high-grade insulation to prevent heat loss. This is provided with ISI mark 1500 watts heating device, controlled by capillary type thermostat from ambient to 110°C, within +1°C., two indicator lights & one On/off Rotary Switch.	
PT-250/05	350 x 250 x 100mm      6 Holes of 75mm	<b>7500-00 each</b>
PT-250/10	400 x 300 x 100mm      12 Holes of 75mm	<b>9950-00 each</b>
PT-250-A	<b>WATER BATH</b> provided with electronic digital temperature controller cum indicator. Size of Inner Chambers    No. of Holes	
PT-250A/05	350 x 250 x 100mm      6 Holes	<b>13950-00 each</b>
PT-250A/10	400 x 300 x 100mm      12 Holes	<b>16950-00 each</b>
PT-250-B	<b>WATER BATH</b> Provided with <b>304 grade Stainless Steel outer body</b> with digital temperature controller.	
	<b>Size of Inner Chambers                      No. of Holes</b>	
PT-250B/05	350x 250x100mm                      6 Holes	<b>17950-00 each</b>
PT-250B/10	400x300x100mm                      12 Holes	<b>19950-00 each</b>
PT-250-C	<b>WATER BATH</b> with original E.G.O. (Imported) Thermostat at an extra cost.	<b>950-00 each</b>
PT-250-D	<i>Spare heater for water bath (ISI mark)</i>	<b>795-00 each</b>

Cat. No.	Item's Description	Rate
PT-255	<p><b>WATER BATH SINGLE WALL</b> Rectangular water bath is suitable for evaporating liquids like alcohol or for melting solids such as agar Concentric rings can be removed to accommodate containers of various sizes.</p> <p><b>CONSTRUCTIONS</b> Single wall chamber made of high quality 304 grade stainless steel provide with ISI mark 1000 watts heating device. This chamber is built in constant level device &amp; having holes of 3" diameter with concentric rings.</p> <p><i>Sizes of chamber                  No. of Holes Dia of Hole</i></p>	
PT-255/05	350x250x100mm                  6                  75mm	4250-00 each
PT-255/10	400x300x100mm                  12                  75mm	5995-00 each
PT-255/15	Spare heater for water bath (ISI mark) (Kettle Type)	950-00 each
PT-257	<p><b>WATER BATH (For Dye Pots)</b> It consists of six stainless steel cups of 250 ml. capacity. It has a thermostatic control, indicator light and a heater of 1.5 KW.</p>	7500-00 each
PT-266	<p><b>CYLINDRICAL OIL BATH</b> Double Walled Provided with Temp. Range 50 to 300°C Controlled by Thermostate</p>	
PT-266/05	S.S. Inner Chamber Cap - 2 Litre	8450-00 each
PT-266/10	Chamber Cap 3/5 Litre	11000-00 each
PT-266/15	Chamber Cap 10 Litre	14500-00 each
<b>OPTIONAL ACCESSORIES</b>		
PT-266/20	Electronic digital temperature controller cum indicator instead of thermostat (extra cost).	7950-00 each
PT-267	<p><b>CYLINDRICAL WATER BATH</b> double wall with stainless steel inner chamber &amp; stand of flask thermostatic controlled.</p>	
PT-267/05	Suitable for 5 litre flask	8500-00 each
PT-267/10	Suitable for 10 litre flask	9950-00 each