

PT-194

TISSUE HOMOGENISERS (Low Speed) It is suitable for low speed homogenisation of tissue and also small volume low speed stirring applications. Low speed operation provides high torque for grinding and homogenising of animal or vegetable tissues. Provided with glass cup with teflon (P.T.F.E.) pestle having stainless steel shaft, capacity 30 ml and variable speed controlled from 200 to 2400 RPM.

12950-00 each



Cat. No.	Item's Description	Rate
PT-196	RUBBER BUSH For stirring (6 or 8 mm.).	120-00 each
PT-200	TEFLON PESTLE TISSUE HOMOGENISERS (GRINDER) Consisting of piston type teflon pestle with serrated tip fitted to a stainless steel rod and a grinding vessel of borosilicate glass for the preparation of tissue homogenates.	
	Working Capacity	
PT-200/05	5 ml.	550-00 each
PT-200/10	10 ml.	650-00 each
PT-200/15	30 ml.	750-00 each
PT-200/20	55 ml.	950-00 each
PT-200/25	100 ml.	1375-00 each

TEFLON MAGNETIC ROTORS

PT-205	TEFLON/P.T.F.E. COATED MAGNETIC ROTORS/NEEDLES	
	Length x Dia.	
PT-205/05	12 x 8 mm.	120-00 each
PT-205/10	18 x 8 mm.	150-00 each
PT-205/15	25 x 8 mm.	175-00 each
PT-205/20	30 x 8 mm.	195-00 each
PT-205/25	35 x 8 mm.	225-00 each
PT-205/30	50 X 9 MM.	275-00 each
PT-205/35	Teflon Policeman 300 mm	650-00 each