



**PT-210
STUFFING BOX GLAND**

Cat. No.	Item's Description	Rate	
TEFLON STUFFING BOX GLANDS			
PT-210	STUFFING BOX GLAND With integral driving pulley, an adequate seal is achieved by adjustable compressed suitable 'O' rings. The cone part is of teflon (P.T.F.E.) and the pulley is made of brass, fitted with ball bearing.		
	<i>Size of Cone</i>	<i>To Fit Stirrer Shaft</i>	
PT-210/05	B-14	6 mm.	1800-00 each
PT-210/10	B-19	6 mm.	1900-00 each
PT-210/15	B-19	8 mm.	1900-00 each
PT-210/20	B-24	6 mm.	2100-00 each
PT-210/25	B-24	8 mm.	2100-00 each
PT-210/30	B-29	8 mm.	2500-00 each
PT-210/35	B-29	10 mm.	2500-00 each
PT-210/40	B-34	8 mm.	2950-00 each
PT-210/45	B-34	10 mm.	2950-00 each
PT-215	STUFFING BOX GLAND All teflon/P.T.F.E. without pulley.		
PT-215/05	B-14	6 mm.	540-00 each
PT-215/10	B-19	6 mm.	575-00 each
PT-215/15	B-19	8 mm.	575-00 each
PT-215/20	B-24	6 mm.	795-00 each
PT-215/25	B-24	8 mm.	795-00 each
PT-215/30	B-29	8 mm.	950-00 each
PT-215/35	B-29	10 mm.	950-00 each

Cat. No.	Item's Description		Rate
	Size of Cone	To Fit Stirrer Shaft	
PT-215/40	B-34	8mm.	1150-00 each
PT-215/45	B-34	10mm.	1150-00 each

TEFLON BLADE STIRRERS

PT-220 **TEFLON (P.T.F.E.) BLADE STIRRERS** The teflon blades can be removed and are mutually interchangeable with the glass shafts, these factors simplify assembly, cleaning and also reduce breakages. These are designed to sweep close to vessel walls and will not cause scratching.

	<i>Approx Shaft Dia. mm.</i>	<i>Shaft Length mm.</i>	<i>Smallest Neck</i>	<i>Blade through Socket</i>	
PT-220/05	6	450	B-19	52	325-00 each
PT-220/10	8	450	B-19	52	340-00 each
PT-220/15	6	450	B-24	76	425-00 each
PT-220/20	8	450	B-24	76	440-00 each
PT-220/25	8	600	B-24	76	465-00 each
PT-220/30	10	450	B-24	76	475-00 each
PT-220/35	8	450	B-29	100	540-00 each
PT-220/40	8	600	B-29	100	565-00 each
PT-220/45	10	600	B-34	100	610-00 each
PT-220/50	10	600	B-34	125	735-00 each
PT-220/55	10	600	B-50	150	935-00 each
PT-221	RUBBER BUSH For Stirring (6 or 8 mm.).				120-00 each
PT-222	TEFLON (P.T.F.E.) BLADED STIRRERS				
	All Teflon Bladed Stirrer 10mm x 450mm with 75 mm Blade.				1100-00 each