



PT - 350 A MUFFLE FURNACE

Cat. No.	Item's Description	Rate										
PT-350	<p>MUFFLE FURNACE Muffle furnaces are suitable for high temperature up to 900°C to fulfill the demands of research institutions, industries & R & D centres.</p> <p>CONSTRUCTIONS : The outer casing is made of thick mild steel finished with epoxy based powder coated paint. Heating elements are made of high temperature kanthal wire surrounding the heating chamber, assuring the most uniform distribution of heat and protects the heating elements against damage and chemical deterioration. Highly efficient lightweight and insulation is cast permanently into the furnace body. All furnaces are equipped with a dependable indicator pyrometer having a 3'5" scale calibrated at 20 °C increments to 1200°C with chromel alumel thermocouple. Thermal fuse, lead wire and plug to work on 220/230 volts AC are also provided.</p> <p>TEMPERATURE : The furnace can be operated continuously up to 900°C temperature controlled by imported sunvic energy regulator fitted in front of the furnace with two indicating lights.</p>											
	<table border="1"> <thead> <tr> <th style="text-align: left;"><i>Muffles Sizes</i></th> <th style="text-align: left;"><i>Ratings</i></th> </tr> </thead> <tbody> <tr> <td>100x100x 225 mm</td> <td>1.6 kw</td> </tr> <tr> <td>125 x 125 x 250 mm</td> <td>2 kw</td> </tr> <tr> <td>150 x 150 x 300 mm</td> <td>3.0 kw</td> </tr> <tr> <td>200 x 200 x 300 mm</td> <td>3.5 kw</td> </tr> </tbody> </table>	<i>Muffles Sizes</i>	<i>Ratings</i>	100x100x 225 mm	1.6 kw	125 x 125 x 250 mm	2 kw	150 x 150 x 300 mm	3.0 kw	200 x 200 x 300 mm	3.5 kw	
<i>Muffles Sizes</i>	<i>Ratings</i>											
100x100x 225 mm	1.6 kw											
125 x 125 x 250 mm	2 kw											
150 x 150 x 300 mm	3.0 kw											
200 x 200 x 300 mm	3.5 kw											
PT-350/05		17950-00 each										
PT-350/10		20950-00- each										
PT-350/15		23950-00 each										
PT-350/20		34950-00 each										
PT-350-A	<p>MUFFLE FURNACE Provided with digital temperature controller cum indicator in lieu of energy regulator and pyrometer.</p>											
	<table border="1"> <thead> <tr> <th style="text-align: left;"><i>Muffle Size</i></th> <th style="text-align: left;"><i>Rating</i></th> </tr> </thead> <tbody> <tr> <td>100 x 100 x 225 mm</td> <td>1.6 k. w.</td> </tr> <tr> <td>125 x 125 x 250 mm</td> <td>2 k. w.</td> </tr> <tr> <td>150 x 150 x 300 mm</td> <td>3.0 k. w.</td> </tr> <tr> <td>200 x 200 x 300 mm</td> <td>3.5 k. w.</td> </tr> </tbody> </table>	<i>Muffle Size</i>	<i>Rating</i>	100 x 100 x 225 mm	1.6 k. w.	125 x 125 x 250 mm	2 k. w.	150 x 150 x 300 mm	3.0 k. w.	200 x 200 x 300 mm	3.5 k. w.	
<i>Muffle Size</i>	<i>Rating</i>											
100 x 100 x 225 mm	1.6 k. w.											
125 x 125 x 250 mm	2 k. w.											
150 x 150 x 300 mm	3.0 k. w.											
200 x 200 x 300 mm	3.5 k. w.											
PT-350A/05		26950-00 each										
PT-350A/10		29950-00 each										
PT-350A/15		33950-00 each										
PT-350A/20		45950-00 each										

PT-350-B **MUFFLE FURNACE** provided with **304 grade Stainless Steel outer body** with digital temperature controller cum indicator.

Muffle Size

Rating

PT-350B/05 100 x 100 x 225 mm 1.6 k. w. **35950-00 each**

PT-350B/10 125 x 125 x 250 mm 2 k. w. **39950-00 each**

PT-350B/15 150 x 150 x 300 mm 3.0 k. w. **43950-00 each**

PT-350B/20 200 x 200 x 300 mm 3.5 k. w. **57950-00 each**

PT-350-C **MICRO PROCESSOR** BASED PID TEMPERATURE CONTROLLER cum indicator in lieu of energy regulator & pyrometer at an extra cost. **12950-00 each**