



PT- 410 A
INCUBATOR BACTERIOLOGICAL

Cat. No.	Item's Description	Rate
PT-410	<p>INCUBATOR BACTERIOLOGICAL (MEMMERT TYPE) Incubators are suitable for variety of uses like drying of slides, araffin embedding, tissue culture work incubation of anti body test, microbiological determination incubation of hydroxy steroids and crystallization studies. Incubators are suitable for temperature up to 80°C with an accuracy of + 1° C.</p>	

CONSTRUCTION All external cases are made of mild steel, often with strong seamless welding & covered with a thick layer of stoved epoxy polyester paint, which provide a strong, highly screech & heat resistant surface for a long lasting attractive finish. Internal surfaces (Including the doors) are manufactured from **304 grade stainless steel**, with great care in order to eliminate sharp edges. The double walled chamber is filled with high - grade insulation to prevent heat loss. Elements are of the sheathed black heat type are put inside beads & placed at the bottom & in both the side ribs are placed for uniform temperature. The door closes on to an asbestos - gasket & fitted with a glass window. The shelves slide neatly into the fixed runners. The fixed vatilators are provided on both side of the wall. The thermometer is provided on the top of incubator to check the temperature. The incubator is mounted on four rubber feet to prevent slipping & this protects the bench surface. The control panel houses a main on/off switch indicator lamp & temperature setting knob.

SIZES

	<i>Width x Depth x Height</i>	<i>Watts</i>	<i>Cap/Liter</i>	<i>Shelves</i>	<i>Stainless steel</i>
PT-410/05	300 x 300 x 300mm.	250	27	1	13500-00 each
PT-410/10	355 x 355 x 355	300	45	2	14500-00 each
PT-410/20	455 x 455 x 455	300	95	2	20950-00 each
PT-410/25	605 x 605 X 605	500	220	3	34950-00 each

PT-410A **INCUBATOR BACTERIOLOGICAL** provided with digital temperature controller cum indicator in lieu of thermostat & fitted with fan to achieve same temperature, in each & every corner of the chamber.

Cat. No.	Item's Description	SIZES				Rate
	<i>Width x Depth x Height</i>	<i>Watts</i>	<i>Cap./Litre</i>	<i>Shelves</i>	<i>Stainless steel</i>	
PT-410A/05	300 x 300 x 300mm.	250	27	1	25500-00 each	
PT-410A/10	355 x 355 x 355	300	45	2	26500-00 each	
PT-410A/20	455 x 455 x 455	300	95	2	33950-00 each	
PT-410A/25	605 x 605 X 605	500	220	3	48950-00 each	
PT-410B	INCUBATOR BACTERIOLOGICAL provided with 304 grade stainless steel outer body with digital temperature controller and air-circulating fan.(G.M.P.)					
	<i>Width x Depth x Height</i>	<i>Watts</i>	<i>Cap./Litre</i>	<i>Shelves</i>	<i>Stainless steel</i>	
PT-410B/05	300 X 300 X300mm.	250	27	1	29950-00 each	
PT-410B/10	355 X 355 X 355	300	45	2	30950-00 each	
PT-410B/20	455 X 455 X 455	300	95	2	43950-00 each	
PT-410B/25	605 X 605 X 605	500	220	3	62950-00 each	
PT-410-C	INCUBATOR BACTERIOLOGICAL provided with P.I.D. Controller instead of Digital Temperature controller cum indicator at an extra cost.					
					5950-00 each	
PT-410-D	INCUBATOR BACTERIOLOGICAL provided with P.I.D. Controller instead of Digital Temperature and RS 232/RS 485 socket for computer attachment you can control the incubator through computer (Desktop).					
					29500-00 each	