



## PT- 6520 SILICONE TUBING

Cat. No.	Item's Description	Rate
PT-6520	<b>SILICONE TUBING (Minimum Packing 5 meter)</b> Silicon tubing is suitable upto -50°C to 250°C this is autoclavable and odourless. So it is suitable for transfer of pharmaceutical liquids beverages and other food preparations. It is not irritant to skin. Due to highly elasticity it is also useful with penestatic pumps.	

	<i>Inner Dia.mm.</i>	<i>Outer Dia. mm.</i>	<i>Wall thickness mm.</i>	<i>Price/Meter</i>
PT-6520/05	2	4	1.0	65-00 each
PT-6520/10	3	5	1.0	95-00 each
PT-6520/15	4	6	1.0	120-00 each
PT-6520/20	5	7	1.0	150-00 each
PT-6520/25	4	7	1.5	150-00 each
PT-6520/30	6	8	1.0	195-00 each
PT-6520/35	6	9	1.5	220-00 each
PT-6520/40	7	10	1.5	275-00 each
PT-6520/45	6	10	2.0	275-00 each
PT-6520/50	8	11	1.5	295-00 each
PT-6520/55	7	11	2.0	295-00 each
PT-6520/60	8	12	2.0	335-00 each
PT-6520/65	10	14	2.0	595-00 each
PT-6520/70	12.5	17	2.5	775-00 each
PT-6520/75	14	20	3	950-00 each
PT-6520/80	16	22	3	1500-00 each