



PT- 6700
LAMINAR AIR FLOW CABINET

Cat. No.	Item's Description	Rate
PT-6700	LAMINAR AIR FLOW CABINET (HORIZONTAL) Supplied complete with stainless steel table top, HEPA filters, pre-filters, U.V. light, static pressure manometer, cock for gas and air or vacuum line.	

	PT-6700/05	PT-6700/10	PT-6700/15	PT-6700/20
Working area	2'x2'x2'	3'x2'x2'	4'x2'x2'	6'x2'x2'
Size of HEPA filter	2'x2'x6"	3'x2'x6"	4'x2'x6"	3'x2'x6"
No. of HEPA filter	1	1	1	2
No. of pre-filter	1	1	1	2
Illumination	1x20W	2x20W	2x40W	2x40W
U.V. Tube	1-1½'x15W	1-1½'x15W	1-3'x30W	2-3'x30W
Prices	48000.00	56000.00	66000.00	85000.00

PT-6710	LAMINAR AIR FLOW CABINET (VERTICAL) Supplied complete with stainless steel table top, HEPA filters, pre-filters, U.V. light, static pressure manometer, cock for gas and air or vacuum line.	
---------	---	--

	PT-6710/05	PT-6710/10	PT-6710/15	PT-6710/20
Working area	2'x2'x2'	3'x2'x2'	4'x2'x2'	6'x2'x2'
Size of HEPA filter	2'x2'x6"	3'x2'x6"	4'x2'x6"	3'x2'x6"
No. of HEPA filter	1	1	1	2
No. of pre-filter	1	1	1	2
Illumination	1x20W	2x20W	2x40W	2x40W
U.V. Tube	1-1½'x15W	1-1½'x15W	1-3'x30W	2-3'x30W
Prices	53000.00	62000.00	72000.00	92000.00

PT-6715 OPTIONAL ACCESSORIES

Exhaust for vertical model

8000-00 each