



PT- 85-86
VACUUM PUMPS

Cat. No.	Item's Description	Rate	
PT-83	FILTER PUMP EDWARD TYPE This mechanism creates vacuum while connecting with water tap (made of extruded brass).	695-00 each	
PT-85	VACUUM PUMP SINGLE STAGE OIL sealed, air cooled, V-belt driven and mounted on base plate with belt guard etc., with motor 220 volts AC single phase.		
	<i>Displacement Litres /Min.</i>	<i>Motor H.P. AC Single Phase</i>	
PT-85/05	32	0.25	14500-00 each
PT-85/10	50	0.25	15950-00 each
PT-85/15	75	0.50	19950-00 each
PT-85/20	100	0.50	22500-00 each
PT-86	VACUUM PUMP DOUBLE STAGE Same as above but with two stages to give better vacuum.		
	<i>Displacement Litres/Min.</i>	<i>Motor H.P. AC Single Phase</i>	
PT-86/05	32	0.25	18500-00 each
PT-86/10	50	0.25	21000-00 each
PT-86/15	75	0.50	24950-00 each
PT-86/20	100	0.50	29800-00 each
PT-86/25	150	1.00	45000-00 each
ACCESSORIES FOR MODEL PT-85 AND PT-86			
PT-87/01	Vacuum gauge with regulator	1600-00 each	
PT-87/05	Vacuum gauge with regulator and non-return valve combined.	2800-00 each	
PT-87/10	Moisture trap with non-return valve and vacuum gauge with regulator combined	6950-00 each	
PT-87/15	Vacustat without mercury range 10mm. to 0.001 mm./Hg. pressure.	5800-00 each	